

## LAB SFE 20ml & LAB SFE 100ml



SFE 100 ml SFE

Two systems are versatile equipment designed for Supercritical Fluid Extraction in the lab.

Particularly adapted for applications where small quantities of products are available, this system can be used to determine the solubility of your products in a supercritical fluid, collect small sample of extracts for analysis and characterisation, screen rapidly a wide range of operating conditions.

For toxic products, this system is small and can be easily contained and placed in a confined environment.

They are designed for an operation up to 500 bar and 150°C but any other specification is available on request.

Optionally, they can be equipped with co-solvent module, particle formation process configuration, gel drying option, sampling loop and analytical equipment, solvent recycling, training by our R&D experts, etc...

### **SEPAREX S.A.**

5 rue Jacques Monod, BP 9, 54250 Champigneulle – France

Tél : +33 (0)3.83.31.24.24. - Fax : +33 (0)3.83.31.24.83.

Société Anonyme au capital de 750 000 € - RCS Nancy B 333 187 649

Siret : 333 187 649 00020 - Code NAF : 742C - Identification TVA : FR 75 333 187 649